

uBioWellness

Your Digital Healthcare Partner in the Palm of Your Hand

Biosense Creative Corporation | www.ubiowellness.com

Contents

- 01 What is uBioWellness?
- 02 Key Functions
- 03 Product Features
- 04 Package Components
- 05 Specifications
- 06 How to Use
- 07 Result Reports
- 08 Application Solutions
- 09 Certifications
- 10 About the Company

| 01 What is uBioWellness?

uBioWellness is
a personal healthcare device designed to
support a healthy lifestyle by analyzing
Biofeedback data.



02 Key Functions



Wellness Check

Analyzes peripheral pulse waves to assess circulation index and vascular flexibility.



Stress Index Check

Analyzes heart rate variability to evaluate stress levels and autonomic nervous system balance.



03 Product Features




Portability



A simple and compact design that fits easily in a handbag or pocket.



Wireless



Synchronization with smartphones and tablets via Bluetooth LE 4.0 (Supports iOS and Android).



Speed



Measures circulation balance and stress index in just 2 minutes and 30 seconds

04 Package Components

Device + Box + User Manual



05 Specifications



Classification	Specification
Model / Manufacturer	uBioWellness / Biosense Creative Corporation
Size / Weight	58 x 33 x 27 (mm) / 36g
Battery / USB-C charging	450mAh Li-ion / USB-C (2 hours charge, 8 hours use)
Wireless Standard	Bluetooth LE 4.0
Supported OS	iOS (iPhone/iPad), Android (Smartphones/Tablets)

06 How to Use

STEP 01 : App Installation



Install the uBioWellness App
via QR code (Android/iOS)

STEP 02 : Bluetooth Connection



Turn on the device and
connect via Bluetooth

STEP 03 : Measurement




Insert finger and select
measurement mode

07 Result Reports


Wellness Report

Visualizes blood circulation and vascular health with scores and grades (A to G).



Stress Index Report

Clear overview of autonomic nervous balance and stress index.



08 Application Solutions



Integrated Wellness Platform based on Big Data and AI

Cloud Data Storage

Secure cloud storage for long-term tracking


AI Customized Analysis

AI-based personalized analysis using big data

Scalability

Flexible integration with mobile, medical, and lifestyle healthcare services

09 Certifications



- > Patents : Domestic patents for pulse wave analysis algorithms
- > Wireless : Bluetooth communication stability
- > KC Certification : EMC & Safety
- > Domestic Production: 100% developed and manufactured in South Korea

| 10 About the Company



Biosense Creative Co., Ltd.

We have over 20 years of experience in IT and BT industries, developing innovative products that enhance quality of life.

 #425, Seoul Technopark, 232 Gongreung-ro, Nowon-gu, Seoul, Korea

 +82-2-960-3170

 biosense@biosense.co.kr

Thank You