




* Transferring data(Mobile Phone or Tablet PC)




My Files -> Internal or Other storage -> Android -> data -> com.biosense.ubiomacpa -> files -> work
-> test report files

(* The following example is the Samsung Galaxy S6.)




* Transferring data(Mobile Phone or Tablet PC)

My Files -> Internal or Other storage -> Android -> data -> com.biosense.ubiomacpa -> files -> work
-> test report files



* Transferring data(Mobile Phone or Tablet PC)

My Files -> Internal or Other storage -> Android -> data -> com.biosense.ubiomacpa -> files -> work
-> test report files



* From Android 11 or higher version, the data is hidden and invisible on mobile phones, so please connect your mobile phone to a PC to back up the data!!